



HAL
open science

**Technical note: Analytical protocols and performance
for apatite and zircon (U-Th)/He analysis on
quadrupole and magnetic sector mass spectrometer
systems between 2007 and 2020**

Cécile Gautheron, Rosella Pinna-Jamme, Alexis Derycke, Floriane Ahadi,
Caroline Sanchez, Frédéric Haurine, Gael Monvoisin, Damien Barbosa,
Guillaume Delpech, Joseph Maltese, et al.

► **To cite this version:**

Cécile Gautheron, Rosella Pinna-Jamme, Alexis Derycke, Floriane Ahadi, Caroline Sanchez, et al..
Technical note: Analytical protocols and performance for apatite and zircon (U-Th)/He analysis on
quadrupole and magnetic sector mass spectrometer systems between 2007 and 2020. *Geochronology*,
2021, 3, pp.351-370. 10.5194/gchron-3-351-2021 . insu-03745270

HAL Id: insu-03745270

<https://insu.hal.science/insu-03745270>

Submitted on 4 Aug 2022

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution 4.0 International License



Technical note: Analytical protocols and performance for apatite and zircon (U–Th) / He analysis on quadrupole and magnetic sector mass spectrometer systems between 2007 and 2020

Cécile Gautheron¹, Rosella Pinna-Jamme¹, Alexis Derycke¹, Floriane Ahadi¹, Caroline Sanchez¹, Frédéric Haurine¹, Gael Monvoisin¹, Damien Barbosa¹, Guillaume Delpech¹, Joseph Maltese², Philippe Sarda¹, and Laurent Tassan-Got²

¹CNRS, GEOPS, Université Paris-Saclay, 91405 Orsay, France

²CNRS/IN2P3, IJCLab, Université Paris-Saclay, 91405 Orsay, France

Correspondence: Cécile Gautheron (cecile.gautheron@universite-paris-saclay.fr)

Received: 17 January 2021 – Discussion started: 25 January 2021

Revised: 19 April 2021 – Accepted: 27 April 2021 – Published: 1 June 2021

Abstract. Apatite and zircon (U–Th) / He thermochronological data are obtained through a combination of crystal selection, He content measurement by crystal heating with analysis using noble gas mass spectrometry, and measurement of U, Th, and Sm contents by crystal dissolution as well as solution analysis using inductively coupled plasma mass spectrometry (ICP-MS). This contribution documents the methods for helium thermochronology used at the GEOPS laboratory, Paris-Saclay University, between 2007 and the present that allow apatite and zircon (U–Th) / He data to be obtained with precision. More specifically, we show that the He content can be determined with precision (at 5 %) and accuracy using a calibration of the He sensitivity based on the Durango apatite, and its use also appears crucial to check for He and U–Th–Sm analytical problems. The Durango apatite used as a standard is therefore a suitable mineral to perform precise He calibration and yields (U–Th) / He ages of 31.1 ± 1.4 Ma with an analytical error of less than 5 % (1σ). The (U–Th) / He ages for the Fish Canyon Tuff zircon standard yield a dispersion of about 9 % (1σ) with a mean age of 27.0 ± 2.6 Ma, which is comparable to other laboratories. For the long-term quality control of the (U–Th) / He data, attention is paid to evaluating the drift of He sensitivity and blanks through time as well as that of (U–Th) / He ages and Th / U ratios (with Sm / Th when possible), all relying on the use of Durango apatite and Fish Canyon Tuff zircon as standards.

1 Introduction

Apatite and zircon (U–Th) / He thermochronology (AHe and ZHe, respectively) is now a mainstream tool to reconstruct the Earth's evolution through the history of cooling and exhumation over the first dozen of kilometers of the Earth's crust (e.g., Farley, 2000, 2002; Gautheron and Zeitler, 2020; Reiners, 2005; Reiners and Brandon, 2006). The geological implications of these data rely on the precision of measurements of He, U, Th, and Sm contents of apatite and zircon crystals by (i) crystal picking and the accuracy of the alpha-ejection correction through crystal shape measurement; (ii) non-destructive He degassing and abundance determination by mass spectrometry; and (iii) U, Th, and Sm analysis after crystal dissolution and solution analysis by inductively coupled mass spectrometry (ICP-MS) as well as on knowledge of He diffusion. Different contributions have already presented parts of the analytical protocols, for example, for crystal dimensions measurements (e.g., Cooperdock et al., 2019; Glotzbach et al., 2019; Herman et al., 2007), He degassing using a laser beam (e.g., Foeken et al., 2006; House et al., 2000), dissolution and analysis of U, Th, Sm, Ca, or Zr (e.g., Evans et al., 2005; Guenther et al., 2016; Reiners and Nicolescu, 2007), and improvements of noble gas analysis by magnetic sector mass spectrometry (Burnard and Farley, 2000). In addition, several contributions have also been dedicated to the determination of He diffusion in apatite and zircon (e.g., Farley, 2000; Flowers et al., 2009; Gautheron et al., 2009, 2020; Gerin et al., 2017; Guenther et al., 2013; Goldsmith et al., 2020; Shuster et al., 2006).

In this contribution, we present all the methodologies that have been developed and used in the GEOPS laboratory, Paris-Saclay University, for more than a decade, focusing on the He degassing process and He content analysis, acid digestion, U and Th (Sm and Ca for apatite) analysis by ICP-MS, and data reduction, which all lead to (U–Th) / He thermochronological data. The detailed protocols represent one method among many others, and this technical note has no other purpose than to present them. Specifically, analytical details are given on two built-in-house noble gas extraction–purification lines, one coupled to a quadrupole mass spectrometer and the other to a magnetic sector mass spectrometer. An efficient method of He content calibration using fragments of the Durango apatite standard is presented. Finally, (U–Th) / He data obtained for 6 to 8 months in the laboratory for Durango apatite and Fish Canyon Tuff (FCT) zircon are presented. The wealth of data, i.e., 272 apatite and 57 zircon measurements, allows us to gain new insights on data acquisition, analytical difficulties, and the reproducibility of He, U, and Th (Sm for apatite) contents as well as (U–Th) / He data.

2 Methods

2.1 Apatite and zircon (U–Th) / He thermochronological methods

The (U–Th) / He thermochronological method is based on the accumulation within the crystal structure of ^4He atoms produced by the decay of ^{238}U , ^{235}U , ^{232}Th , and ^{147}Sm to alpha particles that become ^4He atoms, and it exploits temperature-dependent ^4He diffusion through a crystal lattice. Because the Sm contribution to the He budget is limited, this method is referred to as (U–Th) / He, with (U–Th–Sm) / He used more recently for phases with measurable Sm content (Ault et al., 2019). Specifically, for apatite and zircon, He is well retained within the crystal structure with a temperature sensitivity of ~ 40 – 120 and ~ 20 – 200 °C, respectively (see the recent review of Ault et al., 2019, for details about the complexity of He retention). Upon alpha decay, the alpha particle is emitted with an important kinetic energy, thus traveling along ~ 5 to 30 μm before stopping or being ejected from the crystal (Farley et al., 1996; Ziegler, 2008). Apatite and zircon (U–Th) / He ages thus require an alpha-ejection correction for disintegrations that happen close to crystal margins (Farley et al., 1996). Apatite and zircon (U–Th) / He dating methods then possess specificity in demanding careful extraction of crystals from plutonic, volcanic, or sedimentary rocks using crushing, mineral separation, and crystal selection before He, U, Th, and Sm analyses. Only suitable apatite and zircon crystals should be selected for (U–Th) / He dating; then crystal size and geometry should be measured and recorded in order to correct the alpha loss by ejection.

In this technical note, we are focusing on the protocols developed over more than 10 years at the GEOPS laboratory, Paris-Saclay University, and more specifically on He, U, Th, and Sm content determinations, while some details of the scheme used for apatite and zircon preparation and selection are also given. Contributions that have investigated the impact on (U–Th) / He dating of crystal quality selection, as well as the size and shape of selected crystals, are an important source for those particular points (e.g., Farley, 2002; Brown et al., 2013; Gautheron et al., 2012; Reiners and Farley, 2001; Reiners, 2005).

2.2 Apatite and zircon sample preparation, picking, and packing

Apatite and zircon crystals are extracted from rocks by classical crushing methods, sieved (mesh < 400 μm), and separated following density (tribromomethane and diiodomethane – VWR[®]) and magnetic methods (L-1 Frantz Isodynamic[®] Separator). Mineral separation can also be performed using other less toxic alternatives such as lithium heteropolytungstate, also known as LST. For the Durango apatite gem crystal, a gentle crushing in an agate mortar allows us to obtain fragments of different sizes. Inclusion-free automorphic apatite and zircons or Durango apatite fragments are picked under a binocular microscope (SZX12 – Olympus[®]) and selected as a function of fragment size (> 60 μm). For an automorphic crystal, the length, height, and width are measured, and the termination geometry of the crystal (broken faces, pyramids, no pyramids) is recorded. The ejection factor (F_T) and equivalent sphere radius (R_s , also called ESR) are determined using the Monte Carlo simulation of Gautheron and Tassan-Got (2010), Gautheron et al. (2012), and Ketcham et al. (2011). As the ejection length for the Th decay chain is higher than that for the U decay chain, the value of Th / U ratio is fixed to the measured value to avoid overcorrections or under-corrections (Ketcham et al., 2011; Ziegler, 2008). An internal modification of the Monte Carlo simulation (Qt_LFT software), coupled with an Excel[®] automatic file generation, is used to calculate the F_T , R_s , and crystal weight from a list of different crystal geometries. The Qt_LFT software is available and added as a Supplement to this technical note.

Each apatite and zircon crystal or fragment is respectively placed into a platinum tube (99.95 % purity, 1.0×1.0 mm – Johnson Matthey[®]) or a niobium tube (purity 99.95 %, 1.0×1.0 mm – Alfa Aesar[®]). The choice of the capsule metal (Nb, Pt) for packaging is strictly related to the acid attack protocol. During U, Th, and Sm analysis by means of ICP-MS instruments, the presence of Pt^+ ions at high concentration (> 320 $\mu\text{g mL}^{-1}$) in the sample solution may lead to the formation of the complex platinum argides $^{194}\text{Pt}^{40}\text{Ar}^+$, $^{195}\text{Pt}^{40}\text{Ar}^+$, and $^{198}\text{Pt}^{40}\text{Ar}^+$, which cause isobaric interferences with the measured U isotopes at masses 234, 235, and 238 (Evans et al., 2005; Reiners and Nicolescu, 2007). Ap-

apatite chemical digestion can be achieved using acid digestion at low concentration (HNO_3 5N ultrapure) and low temperature (65°C), which is not able to dissolve the capsule; therefore, Pt tubes can be used. However, an acid digestion using concentrated acids (HCl and pure HF 27N) at high temperature (220°C) is instead used for zircon, which leads to a total dissolution of the Pt capsule. Therefore, Nb tubes were adopted for zircon. Although niobium–argon complexes do not create isobaric interferences with the analyzed U masses, the solution, which is highly concentrated in Nb, may cause a partial precipitation of uranium and thorium (Evans et al., 2005). The Nb impact on the U and Th content determination will be discussed in detail in Sect. 3.3.

2.3 Helium analysis protocol

The helium content analyses were performed at the GEOPS laboratory, Paris-Saclay University (Orsay, France). Each capsule containing a crystal, fragments, or grain(s) was degassed using either a built-in-house He extraction line coupled with a quadrupole mass spectrometer (Prisma QMG 100 Pfeiffer[®]), further referred to as the quad line, or another built-in-house line connected to a rehabilitated VG5400[®] magnetic sector mass spectrometer, further referred to as the VG line. The quad and VG lines are fully automated using LabView[®] software from the heating phase to the helium analysis. Each portion of the line is divided into sections (extraction, purification, analysis) by pneumatic Swagelok[®] valves coupled to electro-valves (E.V 3/2 NF Direct Flasque.D2,4 ALU BUNA, TH France[®]) and activated by pressurized air. Ultrahigh-vacuum conditions (pressure $< 10^{-9}$ mbar) are guaranteed by using a combination of turbomolecular (HighCube – Pfeiffer[®]) and ionic pumps (Star-Cell – Varian[®]). Figure 1 presents the schematic geometry of the two built-in-house quad and VG lines, with their different parts that are controlled using LabView[®].

Platinum and niobium tubes adopted for sample packaging are suitable because of the low level of hydrogen they release under high vacuum, their malleability, and their U, Th, Sm, and rare Earth element (REE) purity requiring no pre-etching. In addition, being metallic materials, they ensure a homogeneous heat transfer during laser shooting. For the quad line, the Pt–Nb tubes are deposited on a copper planchette having 25 or 49 positions, whereas, for the VG line, the tubes are placed on a stainless-steel (inox) or a copper planchette with 12 or 49 positions. Copper or inox was selected for the planchette material due to their good thermal conductivity and their inertia in vacuum conditions. After each sample loading, the line (quad or VG) is heated overnight at low temperature ($< 50^\circ\text{C}$) using heating tapes to remove any gas adsorbed on their inner walls.

For the quad line, the planchette is placed into a cell that moves under the laser beam using an X–Y motor system (SMC100CC – Newport[®]) controlled by LabView[®]. A ytterbium-doped infrared (IR) diode laser (wavelength 1064–

1080 nm, power 10 W, Manlight[®] – Laser2000[®]), coupled with an optic system placed at a focal distance of 4 cm from the sample, allows heating up the capsules with a beam of $\sim 70\ \mu\text{m}$ diameter. For the VG line, the heating of the tubes is ensured by an infrared Fusions Diode Teledyne[®] laser moving above the cell, and ultrahigh vacuum can quickly be obtained by heating an empty capsule purposely placed on the copper–inox planchette with the IR laser to remove air adsorbed on it. In both cases, the cell is sealed with a CF63 sapphire window (Caburn[®]), allowing good transmission of the whole IR laser beam but low He exchange. Each capsule is heated using the heating protocol summarized in Table 1.

The heating schedule procedure is repeated on each sample until all ^4He is degassed, giving a signal back to the background level within less than 2%. The sample temperature achieved using the laser of the quad line is recorded by means of a LabView[®] camera and a built-in-house algorithm that converts the total red, green, and blue visible light into a temperature. We use the visible light emission part of the black-body light emission during heating. To this aim, a Pt capsule was heated with increasing values of laser intensity and pictures were taken at different temperatures, as presented in Fig. 2a. At the same time, the temperature of the heated capsule was measured using an external filament extinction pyrometer. This type of pyrometer is currently used to calibrate thermal ionization mass spectrometer (TIMS) filament temperatures. For each recorded picture, a simple image treatment has been realized using LabView[®] to retrieve the red, blue, and green value on a red–blue–green (RBG) colorimetric coding system that ranges from 0 to 255 (Fig. 2b). As the red signal is already saturated when the capsule is emitting visible light, we chose to sum the signals of the three colors, as shown in Fig. 2c. The obtained RBG sum is correlated with the temperature of the heated capsule. This simple image treatment procedure has been calibrated from 950 to 1150°C for a fixed value of the exposure time of 500 ms. The same image treatment has been automatically applied to each heated capsule, allowing us to retrieve temperature with an estimated error of $\pm 20^\circ\text{C}$ on temperature calculation.

The protocol has been designed for the purpose of degassing apatite crystals or fragments that are packed in a Pt tube; the chosen heating time and temperature schedule should permit us to retrieve all the He from apatite during the first degassing step, as predicted using the mean diffusion coefficient of Farley (2000) for Durango apatite. If the unknown apatite is pure and presents a diffusion coefficient similar to the Durango apatite, the He content of the second step should be similar to background. This protocol thus allows us to detect any He retentive mineral inclusion (e.g., titanite, zircon) that was not seen on picking and contributes to the He budget (Farley, 2002). It can also be used to monitor the impact of radiation damage on the retentiveness of an apatite crystal. For zircon, as mineral inclusions are not an issue for (U–Th) / He age interpretation, we simply heat the Nb capsule at a temperature close to image color satura-




Figure 1. Schematic representation of the built-in-house systems for He extraction, purification, and analysis by a mass spectrometer: (a) quad line and (b) VG line. The part with the He diffusion cell in panel (b) is not used in this study.

tion ($\sim 1150\text{ }^{\circ}\text{C}$; Fig. 2) and avoid any overheating problems such as capsule melting. Because Pt and Nb have different emissivity values, the temperature calibration is not totally adapted to Nb. Indeed, the color visible for a given temperature is not the same for Nb and Pt capsules, and the real temperature of an Nb capsule should slightly differ from that displayed. Nevertheless, the zircon crystals can be heated until total He degassing, and neither the poor determination of

heating temperature nor the different emissivity values are an issue.

The analysis protocols differ as a function of the type of mass spectrometer used and the type of analyzed minerals.

Table 1. Sample packing as well as He purification and analysis protocols.

Minerals	Packing	Degassing	Purification	Mass spectrometer	Standard
Apatite	Pt tube	5 min at 1050 °C*	Liquid N ₂ on activated charcoal + SAES 707 10 min at room temperature	Quadrupole: H ₂ , ³ He, ⁴ He, mass 5, H ₂ O, ⁴⁰ Ar, CO ₂ Electron multiplier 850 V	Volcanic Durango apatite 31.02 ± 1.01 Ma; McDowell et al. (2005)
Zircon	Nb tube	30 min at > 1150 °C*	Liquid N ₂ on activated charcoal + SAES 707 10 min at room temperature	VG5400: ⁴ He Electron multiplier 3500 V + LabView ion counting	Volcanic Fish Canyon Tuff zircon 28.5 ± 0.06 Ma; Schmitz and Bowring (2001)

* Apatite and zircon can be heated at different temperatures and with different time lengths, especially zircon, which should be heated at higher temperature. Equations established by Fechtig and Kalbitzer (1966) can be used to calculate the minimal time for a given temperature that should be used to ensure complete He degassing, knowing the He diffusion coefficient. It is, however, important not to heat apatite crystals at a temperature that is too high to avoid any issues with U or Th volatilization.




Figure 2. Temperature calibration procedure of a heated Pt capsule in visible light. (a) Example of capsules heated with different laser beam intensities showing the change in color. (b, c) Evolution, as a function of the temperature measured using a pyrometer, of the RGB code values and their sum, respectively.

2.3.1 Quad line

The diffused ⁴He gas is mixed in the purification line with a known amount of ³He, which is used as a spike, in concentration of about 100 to 1000 times higher than the ⁴He to be determined. A ~ 4000 cm³ (cubic centimeter) cylinder (V1), filled by ³He gas, is connected to a pipette made by two welded valves with a small, 5 mm diameter, stainless-steel cylinder placed inside to reduce the volume of the pipette (Fig. 1). The approximate volume of the pipette is ~ 0.5 cm³ (V2) and allows an ³He amount of 10⁻⁹ to 10⁻¹⁰ ccSTP (cubic centimeters at standard temperature and pressure) (~ 4 × 10⁻¹⁴ to ~ 4 × 10⁻¹⁵ mol) to be introduced into the line. The amount of ³He decreases in the cylinder by a factor of ~ 0.9999 (V1/(V1 + V2)) for every shot of gas extracted; therefore, the number of pipettes taken is automatically recorded to take the decrease in ³He in the cylinder into account. These statistics allow the data for He age computation to be calculated with the right amount of ³He spike intro-

duced into the line (see also below). According to the spike conditions reported above (³He in concentration of about 100 to 1000 times higher than ⁴He), it has been observed that the quad line can perform a total of about 6000 analyses, after which the ³He gas tank is refilled using a 1 L tank filled with ³He (99.5 % – Eurisotop®) at 1 bar.

The sample gas is purified from most of the H₂O, CO₂, H₂, and Ar gases using two liquid nitrogen-cooled traps of activated charcoal, an St707 SAES® getter operated at room temperature, according to different purification protocols adapted for various minerals (Table 1). The use of a hot (> 850 °C) titanium sponge getter (3–13 mm, 99.95 % metal basis – Alfa Aesar®) is dedicated to minerals with high CO₂ or H₂O contents. The access to the entire system of traps, individually connected to the line by means of ultrahigh-vacuum valves (Fig. 1), allows the analysis of a large variety of minerals containing variable abundances of CO₂ or H₂O, such as calcite or goethite (Allard et al., 2018; Cros et al., 2014). Besides helium isotopes (³He and ⁴He), H₂O,

CO₂, H₂, and Ar gases are additionally measured on the electron multiplier of the Prisma QMG 100 Pfeiffer® quadrupole mass spectrometer to check for the effective purification of the analyzed gas. Such measurements of the gas are repeated 16 times, and a linear regression of the data for the ⁴He / ³He ratio is then calculated and includes a correction of the HD⁺ isobaric contribution to the ³He signal, even if this contribution is insignificant compared to the ³He spike signal. In addition, we also observe that the signal at mass 4 slightly increases when the H₂ signal is higher, which we interpret as either a double H₂ molecule having an isobaric impact on mass 4 or the tail of the H₂ peak having an influence on the shape of the mass 4 peak. This effect is not negligible for the low ⁴He signals of typical samples, but an adapted H₂ purification protocol allows us to remove this effect: an St701 SAES® getter unit is positioned just ahead of the quadrupole mass spectrometer (Fig. 1). A longer gettering time could be offered as a suitable alternative to reduce the H₂ influence on the ⁴He peak if the addition of a getter to the mass spectrometer is impossible.

The gas purification protocol, combined with ³He spiking, ensures a close-to-constant total pressure in the line and in the quadrupole mass spectrometer. The ⁴He abundance is calculated from the introduced ³He amount as

$${}^4\text{He} = \left(\left(\frac{{}^4\text{He}}{{}^3\text{He}} \right)_s - \left(\frac{{}^4\text{He}}{{}^3\text{He}} \right)_b \right) \times {}^3\text{He}, \quad (1)$$

with (⁴He / ³He)_s and (⁴He / ³He)_b being the isotopic ratios measured for the sample and blank, respectively. The ³He amount is determined from the calibrated value, denoted ³He_c, the two volume values (*V*1 and *V*2), and two evolution parameters (*N* and *D*) using the following equation:

$${}^3\text{He} = {}^3\text{He}_c \times \left(\frac{V1}{V1 + V2} \times D \right)^N, \quad (2)$$

where *V*1 is the volume of the ³He cylinder (~4000 cm³), *V*2 the pipette volume (~0.5 cm³), *N* the pipette number (i.e., *N* is the number of introductions of the pipette volume), and *D* the “drift”, an additional parameter introduced to account for the evolution of the sensitivity of the quadrupole mass spectrometer along with external parameters such as temperature or even power failures. Both parameters ³He_c and *D* are obtained empirically by calibration from several analyses of Durango apatite fragments. *D* acts as if the pipette volume *V*2 could vary to mimic the variations of the quadrupole sensitivity and is updated manually according to the results of the Durango standard, which is analyzed regularly, especially every time the source is tuned again. The product ³He_c × ($\frac{V1}{V1+V2} \times D$)^{*N*} thus decreases regularly and homogeneously.

2.3.2 VG line

The diffused ⁴He is purified from the H₂O, CO₂, H₂, and Ar gases using two St707 SAES® getters and a Ti sponge getter (Fig. 1b). A cryogenic trap from Advanced Research Systems® (ARS), installed more recently, has the capacity to cool activated charcoal down to 8 K. At this temperature, He is efficiently trapped, then further released into a smaller volume at about 50 K. Again, according to the nature of the sample to be analyzed, different purification protocols are adopted (Table 1). Every protocol is fully automatized, and the ⁴He gas is introduced into the VG5400® magnetic sector mass spectrometer. The filament amperage is chosen to obtain a compromising value for the trap current of 300 to 400 μamps, which is lower than the recommended value for He analysis but ensures a longer lifetime of the filament. Isotope ⁴He is analyzed by ion counting using a Pfeiffer® electron multiplier (17 dynodes) connected to an Ortec® discriminator and a LabView® counting card. A total of 20 analyses of the ⁴He signal integrated over 1 s are performed, and the mean ⁴He signal is recorded (while a linear regression is available, only the mean signal and associated standard deviation are recorded). The dead time of the electronic chain is close to the width of the pulse delivered by the electron multiplier, which is a few nanoseconds. The maximal recorded counting rate being about 3 × 10⁵ cps (counts per second), the dead time correction is always lower than 1 % and is neglected. The system sensitivity is empirically determined from ⁴He analyses of the Durango standard (see Sect. 3.2), while the internal ⁴He standard from a ~4000 cm³ cylinder is only used to check for sensitivity change over time and signal reproducibility. The ⁴He abundance in a sample is computed using

$${}^4\text{He} = \left({}^4\text{He}_s - {}^4\text{He}_b \right) \times S, \quad (3)$$

where ⁴He_s and ⁴He_b are the signals measured for the sample and blank, respectively, and *S* is sensitivity in ccSTP He cps⁻¹. The unit ccSTP is not the SI unit for the amount of a substance, but it is a way of expressing gas moles by the volume (cubic centimeter) they would occupy at the standard temperature of 273.15 K = 0 °C and the standard pressure of 101 325 Pa = 1 atm. This definition does not follow the 1 bar pressure 1982 IUPAC recommendation but is historically used by the noble gas community. We thus keep this definition of standard pressure and temperature conditions, STP, with 1 mol occupying 22 414 ccSTP.

Durango apatite fragments and/or Fish Canyon Tuff zircon crystals are analyzed regularly (one Durango or Fish Canyon standard analyzed every seven unknown samples) to check the (U–Th) / He analysis reproducibility.

2.4 Digestion chemistry protocol

2.4.1 Vessel cleaning

For apatite, we used single-use 4 mL polypropylene (PP) snap-cap tubes (supplier VWR[®]) that do not need prewashing for our purpose (Table 2).

For zircon, the dissolution is made in 350 μL PFA (perfluoroalkoxy Teflon[®]) Parrish-style vials (Savillex[®]) placed into a high-pressure–high-temperature dual-wall digestion vessel. Before their use for acid digestion of samples, the vials undergo a series of acid baths in a 250 mL beaker (borosilicate glass – VWR[®]) placed on a hot plate at 100 °C, according to the following sequence. Cycle 1: 24 h bath in diluted 5 % Extran[®] MA 02 (Merck[®]) in Milli-Q[®] water (Milli-Q[®] HX 7000 SD); cycle 2: 24 h bath in HCl 5N (Emsure[®] 32 % VWR[®]); cycle 3: 24 h bath in HCl/HNO₃ 3 / 1 (HCl Emsure[®] 32 % VWR[®]; HNO₃ Emsure[®] for analysis 65 % VWR[®]); cycle 4: 24 h bath in HNO₃ 5N (Emsure[®] 65 % VWR[®]). Between each bath, vials are rinsed with Milli-Q[®] water (18 M Ω Direct 8 System – Merck[®] Millipore[®]). Vials are finally dried in an oven at 50 °C and stored in a hermetically closed PP box until further use (Table 2).

To check their degree of cleanness after the series of cleaning baths, cleaned vials are filled with 200 μL Milli-Q[®] water and refluxed for 2 h at 100 °C; this water is transferred to PP tubes with the addition of 800 μL Milli-Q[®] water for analysis by ICP-MS of ²³⁸U and ²³²Th natural isotopes as well as spike isotopes ²³⁰Th and ²³⁵U.

2.4.2 Spike solution composition and calibration

According to the isotopic dilution method, for every apatite sample to be analyzed, a volume of 50 μL of spike (solution MR2) is introduced into the vial before the dissolution protocol (100 μL for zircon). This spike solution is prepared (60 mL) every 6 to 12 months from elemental mother solutions MR1 and MR (Appendix A, Table A1), these being obtained from certified concentrated mono-elemental solutions: for MR, we mix ²³⁵U (10 mL, IRRM-50, 4.2543(11) nmol ²³⁵U g⁻¹, EU JRC – Belgium), ²³⁰Th (5 mL, IRRM-61, 2.474(18) nmol ²³⁰Th g⁻¹, EU JRC – Belgium), and ¹⁴⁹Sm (100 mL, 7 ppm, 90 %, 47.01 nmol ¹⁴⁹Sm g⁻¹, Berry and Assoc[®] – MI, USA). For apatite analysis, ⁴²Ca (2.44 mg solid ⁴²CaCO₃, 97.8 %, BuyIsotope[®] – Sweden) is added. Detention of EU JRC materials is authorized under laboratory license for radioactive materials (RMs). Other isotopes, such as ²³³U or ²²⁹Th, are also available as isotopic spikes at JRC; however, their high cost due to radioactive nuclear material transport fees is not negligible. For this reason, we chose to use isotopes ²³⁵U and ²³⁰Th. In more than 10 years of research at the GEOPS laboratory, 16 solutions of spike for apatite and three solutions of spike for zircon have been produced. Since the early development of the (U–Th) / He dating method at GEOPS, isotopes other than ²³⁵U and ²³⁰Th have recently been added to the

protocol, namely ¹⁴⁹Sm and ⁴²Ca. For zircon, the addition of the isotope ¹⁴⁹Sm to the spike solution was discarded due to the negligible impact of ¹⁴⁷Sm decay on the He budget compared to the He production from ²³⁸U and ²³²Th.

Before use, the spike solution MR2 is calibrated for U, Th, Sm, and Ca using a series of weighted solutions, each of them obtained by mixing a volume of another mono-elemental standard solution with a volume of the spike MR2 (Appendix A, Table A1): the U_{SS} solution contains a ²³⁸U mono-elemental solution at 1015 ppm (99.96 %, Analab[®] France), diluted and mixed with the ²³⁵U spike; the Th_{SS} solution contains a ²³²Th mono-elemental solution at 993 ppm (99.93 %, Analab[®] France), diluted and mixed with the ²³⁰Th spike; the Sm_{SS} solution contains a natural Sm solution at 1006 ppm (¹⁴⁷Sm = 14.99 %, ¹⁴⁹Sm = 13.82 %, Analab[®] France), diluted and mixed with the ¹⁴⁹Sm spike; and the Ca_{SS} solution contains a natural Ca solution at 1003 ppm (⁴²Ca = 0.65 %, ⁴³Ca = 0.13 %, Analab[®] France) mixed with the ⁴²Ca spike. These solutions are then analyzed by ICP-MS to properly calibrate the MR2 spike solution.

2.4.3 Sample digestion protocol

After degassing, each Pt–Nb-conditioned sample is transferred from the planchette into a vial for grain dissolution by acid digestion. Sample digestion is started by adding a volume of 50 μL for apatite (or 100 μL for zircon) of a spike solution into each vial (in HNO₃ 5N and containing a known amount of ²³⁵U and ²³⁰Th, plus addition of ¹⁴⁹Sm and ⁴²Ca for apatite). According to the nature of the sample, a specific dissolution protocol is followed (acids and heating temperature) (Table 2).

For apatite, the dissolution protocol (Table 2) was adapted from Farley (2002). Dissolution requires a soft acid digestion (HNO₃ 5N – bidistilled from HNO₃ 65 – Normapur[®] – VWR[®]) performed in a 4 mL single-use polypropylene tube (VWR[®]) by adding 50 μL of spike (\sim 4 ppb of ²³⁵U, ²³⁰Th, ¹⁴⁹Sm, and \sim 1 ppm ⁴²Ca) and 50 μL of HNO₃ 5N. The tube is then placed on a hot plate at 65 °C during 3 h for digestion. After digestion and cooling, samples are diluted with 1.9 mL of HNO₃ 1N and stored at 4 °C before ICP-MS analysis (Table 2). Due to the digestion conditions (using diluted acids), the Pt capsule does not dissolve and does not interfere with ICP-MS analysis. Sample digestion is always made with nitric acid, freshly diluted by addition of Milli-Q[®] water. To minimize possible contamination sources from the environment, storing tubes, or evaporation of the solutions, ICP-MS analysis is scheduled within a few days of sample digestion. After analysis, the Pt capsules are promptly collected, cleaned, and sent back to the manufacturing company according to a recycling or reselling loop contract.

For zircon, the dissolution protocol was slightly adapted from Reiners (2005) and Reiners and Nicolescu (2007). The dissolution is performed in 350 μL PFA Parrish-style vials (Savillex[®]). Zircons are first spiked with 100 μL of ²³⁵U and

Table 2. Cleaning chemistry and chemistry protocols for prepared solution.

Mineral	Vial type	Cleaning protocol	Spikes	Dissolution protocol	Solutions
Apatite	4 mL single-use PP tube	No cleaning needed	50 μL ($^{235}\text{U} \sim 4$ ppb; $^{230}\text{Th} \sim 4$ ppb, $^{149}\text{Sm} \sim 4$ ppb; $^{42}\text{Ca} \sim 1$ ppm)	+ 50 μL 5N HNO_3 and heating for 3 h at 65 °C. + cooling time (30 min) + 1.9 mL 1N HNO_3	Blk: acid blank Blk-ch: acid blank chemistry BSP: spiked blank BSP-ch: spiked blank chemistry
Zircon	350 μL PFA tubes + PFA container + stainless-steel digestion bomb + oven + 4 mL single-use PP tube	PFA vials on a hot plate at 100 °C: Extran 5 % bath (24 h); HCl 5N Emsure bath (24 h); HCl: HNO_3 (1 : 3) Emsure bath (24 h); HNO_3 5N Emsure bath (24 h)	100 μL ($^{235}\text{U} \sim 45\text{--}55$ ppb; $^{230}\text{Th} \sim 15\text{--}20$ ppb)	Step 1: vial: + 200 μL HF 27N + few drops HNO_3 7N vessel: 10 mL HF 27N + 1 mL HNO_3 7N oven: 220 °C for 96 h Step 2: complete evaporation at 100 °C Step 3: vial: +300 μL HCl 6N; vessel: 12 mL HCl 6N; oven: 220 °C for 24 h Step 4: complete evaporation at 100 °C Step 5: add 200 μL HNO_3 5N + some drops HF 0.1N; 1 h reflux at 100 °C Step 6: transfer to PP vial, add 800 μL HNO_3 1N Step 7: transfer to PP vial, dilution 1/10 with HNO_3 1N	BSP-ch-Pt: spiked blank chemistry with a Pt capsule for apatite BSP-ch-Nb: spiked blank chemistry with an Nb capsule for zircon DUR: natural Durango Sp: spiked sample including standard

^{230}Th (~ 45 to ~ 55 ppb of ^{235}U , ~ 15 to ~ 20 ppb of ^{230}Th ; Table 2). The vials are then placed into a digestion vessel hermetically sealed with a metallic gasket to hold high pressures (IN/PFA-OUT/stainless-steel PA4748, Parr Instrument Company[®]). The digestion follows several steps summarized in Table 2. Step 1: inside the vials, add 200 μL HF 27N (Suprapur[®] – VWR[®]) and a few drops of HNO_3 7N (Suprapur[®] – VWR[®]); inside the digestion vessel, add 10 mL HF 27N and 1 mL of HNO_3 7N. Once sealed, the vessel is heated up at 220 °C in an oven and held at high pressure for 96 h. Step 2: the acid solution is evaporated to dryness by placing the vials on a hot plate at 100 °C. Step 3: 300 μL HCl 6N (Suprapur[®] – VWR[®]) is added to each vial; the Parr vessel is filled with 12 mL HCl 6N, sealed, and heated back at 220 °C under pressure for 24 h in an oven. Step 4: the vials are evaporated to dryness on a hot plate at 100 °C. Step 5: a reflux is carried out with a combination of HNO_3 5N (200 μL) and HF 0.1 N (a few drops) at 100 °C for 1 h. Step 6: the solutions are transferred to 4 mL polypropylene tubes where 800 μL HNO_3 1N is added. Step 7: a final dilution (1/10) is done with freshly prepared HNO_3 1N in a second 4 mL PP single-use tube. The solutions are then stored at 4 °C before isotopic analysis. To avoid pollution released from the storage tubes or changes in concentration by evaporation of the solutions, the ICP-MS session is always scheduled within a few days after sample dissolution.

In addition to the apatite and zircon samples, the following types of solution are also prepared, as summarized in Table 2: spiked sample (Sp), including Durango and FCT standards; acid blank (Blk) to check acid purity and potential contaminations of tubes; chemistry acid blank (Blk-ch) to check the enrichment contamination in acid caused by the chemistry protocol; spiked blank (BSP) with a weighted vol-

ume of spike added to a volume of acid in order to check the variations in concentration of the spike and to take the contribution of natural isotopes contained in the spike into account; spiked blank chemistry (BSP-ch) with a volume of spike in acid undergoing the same dissolution protocol than the samples, which allows us to quantify contamination coming from contingencies during the dissolution protocol (vessel, user, acid); Durango solution (DUR) with a single fragment of Durango dissolved in a volume of acid and no spike added, which allows the natural isotopic ratio of uranium to be measured for checking, as it has to be in isotopic equilibrium with $^{238}\text{U} / ^{235}\text{U} = 137.88$ (or $^{235}\text{U} / ^{238}\text{U} = 0.00725268$).

2.5 U and Th (Sm and Ca) analysis by ICP-MS

The solutions obtained by chemical dissolution of the samples are then analyzed by ICP-MS to determine the U, Th, and Sm (Ca for apatite) signal intensities. Since 2016, we have mainly used a high-resolution inductively coupled plasma mass spectrometer (HR-ICP-MS – ELEMENT XR[®] from Thermo Scientific[®]) at the GEOPS laboratory that allows the U, Th, and Sm isotopes to be measured at low resolution (300), while Ca is measured at high resolution (10 000). In addition, we also use a quadrupole ICP-QMS series II CCT Thermo-Electron[®] at LSCE (Gif/Yvette; France) and an Agilent[®] 7900 quadrupole ICP-MS at IPGP (France) to measure the U, Th, and Sm contents.

The U, Th, Sm, and Ca abundances are then deduced from the measured $^{235}\text{U} / ^{238}\text{U}$, $^{230}\text{Th} / ^{232}\text{Th}$, $^{149}\text{Sm} / ^{147}\text{Sm}$, and $^{42}\text{Ca} / ^{43}\text{Ca}$ isotopic ratios and equations derived from Evans et al. (2005) by means of a built-in-house Excel[®] workbook, which is available upon request to Alexis Derycke, based on the VBA automation software. Appendix B presents the equations used for the U, Th, Sm, and Ca content determi-

nation. The measurement of the Ca content and the deduced apatite crystal weight, combined with the measurements of U, Th, and Sm abundances, then allow the U, Th, and Sm concentrations of the Durango fragments to be determined. Determination of crystal weight is also useful to ensure that the criteria imposed for grain selection, i.e., (i) crystal size ($L, W, T > 60 \mu\text{m}$), (ii) geometry (well-shaped prisms), and (iii) purity (inclusion-free crystals), have been respected.

2.6 (U–Th) / He age reduction

For a sample, the (U–Th) / He age (Ma) is calculated assuming a linear production of ^4He with time using the determined U, Th, and Sm abundances and the equations presented in Appendix C. The assumption that the ^4He production is linear holds for ages lower than ~ 150 Ma, since the half-lives of the U, Th, and Sm isotopes are ~ 5 to ~ 600 times longer than 150 Ma. When the age is older than 150 Ma, we use a trial-and-error code to calculate the age, which never changes by more than a few percent.

For an ICP-MS session, the (U–Th) / He data reduction Excel[®] workbook can calculate the U–Th–Sm concentrations (ppm by mass) and their associated standard deviations for the different minerals analyzed, according to the chemical dissolution protocol followed, as well as the effective uranium content, eU (ppm), the Th / U and Sm / Th ratios, and the (U–Th) / He ages.

3 (U–Th) / He results and discussion

3.1 Helium quadrupole analysis

One aspect of quadrupole mass spectrometry is the variable response in terms of the ionization and signal generation of such instruments. This behavior has also been observed on the quadrupole adopted here for analyses of rare gases and is illustrated in Fig. 3a. The signal of the ^3He spike is reported as a function of the number of extracted pipettes over a period of 5 months of analysis. The ^3He signal fluctuates significantly although the ^3He amount in each pipette should decrease smoothly following a law that depends on the volume of the cylinder and the pipette volume (here ~ 4000 and $\sim 0.5 \text{ cm}^3$, respectively; see Eq. 2). The advantage of using a ^3He spike for isotopic dilution is to thwart the impact of the nonlinear answer of the quadrupole mass spectrometer with time (Farley, 2002; House et al., 2000), but it also allows the total gas pressure to be buffered in the mass spectrometer if the introduced ^3He signal is large enough compared to the other signals.

Figure 3b presents the various signals measured with the quadrupole mass spectrometer (H_2 , ^3He , ^4He , ^{40}Ar , CO_2 , and mass 5 that represents the background noise) during the same 5 months of analyses (~ 10 to 30 analyses per day, 5 d a week). The ^3He clearly controls the total pressure inside the mass spectrometer during the analysis independently from

the pressure of ^4He gas released from the samples. The use of a spike rich in ^3He allows us to achieve a stable and uniform total pressure in the mass spectrometer for any degassed sample analyzed. In addition, we follow the H_2 , CO_2 , and Ar signals to guard against any leak or gas purification problem. For example, in Fig. 3b, the H_2 signal was higher than usual at a pipette number of 200 to 250, alerting us about a technical problem related to purification. It will then help to take the AHe ages with caution for low ^4He samples, considering the H_2 tail can influence the ^4He signal (Sect. 2.3.1).

However, the quadrupole signals are variable and sometimes erratic within a few percent of the signal over weeks to months (Figs. 3a and 4a). To correct for this quadrupole drift, we introduce a correcting factor, denoted D , and a calibrated ^3He value, denoted $^3\text{He}_c$, in the calculation of the ^3He pipette content (Eq. 2), both of which we adjust using the deviation of the Durango ages obtained over a period of 1–2 months. Parameter D can be viewed as modifying the value of the pipette volume, V_2 in Eq. (2), to counterbalance the variations of the quadrupole sensitivity. As we determine the U and Th (Sm and Ca for apatite) contents every 1 to 2 months, we also determine the D and $^3\text{He}_c$ values for each batch of analyses (unknown samples and Durango standards) when we are calculating the (U–Th) / He ages.

Figure 4c and d present the Durango apatite (U–Th) / He ages obtained from measurements over a period of 2 months using two coupled values of the drift parameter, D , and the $^3\text{He}_c$ calibrated content. We choose to present this example to explain how we perform the calibration and because of the presence of age outliers. Firstly, one can observe that the (U–Th) / He ages obtained for the Durango apatite remain constant and are more reproducible using a D value of 99.937×10^{-2} compared to 99.987×10^{-2} . In the present case, it just depletes the tank in ^3He more quickly. We consider, over the analysis batch time (here pipette number 300 to 525), the variation of the ^3He signal and the lowest standard deviation of Durango AHe ages to empirically determine the D and $^3\text{He}_c$ values to be used (Fig. 4a, b, c, and d). Each pair of values of the $^3\text{He}_c$ and D parameters are recorded, allowing us to detect any problems during the measurements. Most of the time, the sensitivity does not evolve and the same D and $^3\text{He}_c$ values are used. However, each time the filament or the voltage of the multiplier is changed, those values need to be tuned. This protocol allows us to achieve better reproducibility for ages of Durango and unknown samples by 1 % to 3 %. This is not significant in comparison to the dispersion of AHe and ZHe ages of most natural samples but allows us to obtain more reproducible ages. In addition, this way of determining the mass spectrometer sensitivity makes precise knowledge of the volumes of the pipette system (V_1 and V_2) not so important.

Secondly, we also check for any analytical issues with the apatite (U–Th) / He ages, including Durango apatite, during the analysis batch. As observed in Fig. 4c–d, three Durango ages are lower or higher than the expected age of




Figure 3. Evolution of the measured ^3He signal with pipette number (from January to June 2019). **(a)** Decrease in the spike ^3He signal. **(b)** Evolution of all the measured signals (^3He : black dots; ^4He : open black dots; H_2 : orange dots; mass 5 (background noise): light blue dots; ^{40}Ar : green dots; CO_2 : dark blue dots; crosses: summation of all the signals).




Figure 4. Evolution of the Durango apatite age as a function of pipette number. **(a)** Evolution of the ^3He signal over the period of analysis. **(b)** Evolution of the calculated ^3He content of each pipette for two different D drift values (black dots or red dots) and $^3\text{He}_c$ calibrated content using Eq. (2). **(c, d)** Calculated Durango AHe ages for the different D drift values and $^3\text{He}_c$ calibrated content.

31.02 ± 1.01 Ma from McDowell et al. (2005). Those outliers can indicate that the Pt tube was too closed, preventing the acid from penetrating the closed tube, or it can indicate an issue during acid digestion, an analytical issue during ICP-MS analysis, a heating temperature that is too high and has vaporized U and Th, or an issue during He quadrupole analysis. Because, in quadrupole gas analysis, we record the tube temperature during heating of the sample and backgrounds, the issue should be related to a digestion problem or the ICP-MS analysis. For the cases presented in Fig. 4c–d, the lower ages were associated with Th / U ratios that are significantly different from other aliquots, tracing the incomplete recovery of U, Th, and Sm associated with non-total digestion or clog-

ging of the solution during ICP-MS analysis. We thus check unknown samples for any digestion or ICP-MS issues.

In addition to the ^3He pipette number associated with each measurement, we use a specific code name. As an example, for a Durango apatite, the code name D19P11A can be read as Durango, year 2019, planchette no. 11, aliquot A. This designation using the name of the sample and its analysis year allows us to better organize the He analysis database and data backup from a chronological point of view.

3.2 Helium magnetic sector analysis

In comparison to quadrupole mass spectrometers, magnetic sector mass spectrometers have a more stable and linear re-

sponse of ionization and thus allow for better analysis. To test the response of the modified VG5400[®] mass spectrometer, determine the ⁴He sensitivity, and calibrate the ⁴He cylinder, we perform multiple analyses of fragments of Durango apatite having different sizes (dozens to hundreds of micrometers). After degassing, the U, Th, and Sm contents (and AHe ages) of the fragments are determined, and, assuming an age of 31.02 ± 1.01 (McDowell et al., 2005), the amount of ⁴He is calculated for each fragment. The calculated ⁴He contents are compared to the measured ⁴He signals in Fig. 5a, showing a correlation between the calculated ⁴He in ccSTP and the ⁴He signal in counts per second (cps), allowing the sensitivity to be determined. The ratio of the calculated ⁴He to the measured ⁴He normalized to the obtained sensitivity is plotted versus measured ⁴He in Fig. 5b. This operation was performed in 2016 and 2018, as the VG5400[®] was tuned with different source parameters and using different line conditions: since 2018, we have used a cryogenic trap, allowing us to concentrate the gas in a smaller volume and then increase sensitivity (Fig. 5a).

For the two different conditions, the measured ⁴He signals and the ⁴He contents calculated from measured U–Th–Sm display a very good linear correlation ($r^2 > 0.99$, Fig. 5a). Derived sensitivity values, S , were 1.1×10^{-12} ccSTP He cps⁻¹ (4.9×10^{-17} mol He cps⁻¹) in 2016 and 2.9×10^{-14} ccSTP He cps⁻¹ (1.3×10^{-18} mol He cps⁻¹) in 2018, additionally showing that the use of the cryogenic trap increases the sensitivity by a factor of ~ 40 . Figure 5b also reveals the dispersion of the calculated ⁴He content. For the 2016 dataset, a large range of Durango apatite fragments were selected, including smaller fragments of $< 50 \mu\text{m}$ in size, to get a large range of ⁴He signals. However, this may introduce difficulty associated with U, Th, and Sm content determination in very small grains. The dispersion of the normalized ratio of calculated content to the He signal observed for a low ⁴He signal (Fig. 5b) may reflect the difficulty of determining the U–Th content for very small Durango fragments. For the 2018 sensitivity calibration, we selected larger Durango fragments that reduce data dispersion (Fig. 5b).

Such analyses of Durango apatite fragments are also useful to test the electron multiplier and counting system responses, which turned out to be linear from thousands to hundreds of thousands of cps for He, without any impact of the dead counting time. However, we always try to select standard fragments yielding signals similar to unknown samples. Further, the use of both a ⁴He cylinder and Durango apatite fragments is important to follow the evolution of the filament and analyzer conditions through time. The ⁴He tank thus gets calibrated and is mainly used to check for signal reproducibility change and analyzer stability. In addition, Durango apatite fragments are analyzed regularly to follow the sensitivity evolution through time, particularly after each power failure, that cause slight modifications of the filament emission condition.

3.3 U, Th, and Sm chemistry and blanks

Acid blanks are regularly analyzed allowing the acid quality to be controlled. Low intensities are measured for ²³⁵U and ²³⁰Th, i.e., < 20 cps. Meanwhile, the intensities for ²³⁸U and ²³²Th are 100 times higher. Higher signal intensities for all the isotopes are measured in Blk-ch, indicating that the chemical dissolution method adds some contamination to the sample solutions. Such a contribution is nevertheless very negligible compared to the intensities of the signals observed for the apatite and zircon samples (100 000 cps).

For apatite, U, Th, and Sm in blanks are low in comparison with the U, Th, and Sm contents of the apatite, as already stated by Reiners and Nicolescu (2007). Figure 6 presents the evolution of the measured ²³⁵U / ²³⁸U, ²³⁰Th / ²³²Th, and ¹⁴⁹Sm / ¹⁴⁷Sm isotopic ratios for spiked blank (BSP), spiked blank chemistry (BSP-ch), and Durango apatite sample (Sp) solutions. For a spiked sample (Sp), the ²³⁵U / ²³⁸U, ²³⁰Th / ²³²Th, and ¹⁴⁹Sm / ¹⁴⁷Sm ratios range between the BSP-ch value and the natural value (i.e., ²³⁵U / ²³⁸U = 0.00725268; ²³⁰Th / ²³²Th = 0, with no natural ²³⁰Th atoms, and ¹⁴⁹Sm / ¹⁴⁷Sm = 0.9), as can be observed in Fig. 6. The isotopic ratio values for the BSP and BSP-ch blanks do not vary by more than a few percent through the different analyses and are orders of magnitude higher than for the sample (Sp). The BSP-ch displays lower ²³⁵U / ²³⁸U, ²³⁰Th / ²³²Th, and ¹⁴⁹Sm / ¹⁴⁷Sm values compared to the BSP (Fig. 6), showing that the chemistry protocol has an impact on the U, Th, and Sm isotopes in solution. However, this effect remains insignificant compared to the Durango apatite U, Th, and Sm contents and does not influence the (U–Th) / He ages. Nevertheless, we always take care to have well-characterized blanks because for natural apatite crystals the U, Th, and Sm contents are usually lower than for Durango.

In opposition to apatite, the chemical dissolution for zircon requires the acquisition and maintenance of more complex laboratory material, such as PFA vials, a high-temperature, high-pressure Parr[®] bomb, and concentrated acids like HF. Tests on all types of blanks (BSP, BSP-ch, and BSPch-Nb) have been carried out over 18 series of sample analyses performed between March 2016 and February 2018. Under constant vessel cleaning protocol and acid quality conditions, comparison of the analyzed U and Th isotopes is presented in Fig. 7. Isotopic ratios measured in BSP-ch (²³⁵U / ²³⁸U = 16.1 ± 0.3 ; ²³⁰Th / ²³²Th = 4.8 ± 0.5) showed a slight reduction of their values compared to those found in BSP (²³⁵U / ²³⁸U = 17.3 ± 1.0 ; ²³⁰Th / ²³²Th = 5.1 ± 0.2), again indicating a small contamination during the chemical protocol (Fig. 7). The loss of elements and more specifically of Th when HF+HNO₃ acids are used has already been well described (e.g., Révillon and Hureau-Mazaudier, 2009; Yokoyama et al., 1999), and the effect can be seen in Fig. 7, where the ²³⁵U / ²³⁸U and ²³⁰Th / ²³²Th ratios of BSP-ch are lower than for




Figure 5. VG5400[®] magnetic sector mass spectrometer sensitivity determined using fragments of Durango apatite. **(a)** Log–log diagram presenting the calculated ^4He content (nccSTP) versus the measured ^4He signal (cps: counts per second) for two determinations in 2016 and 2018. **(b)** Ratio of the calculated to measured ^4He normalized to the obtained sensitivity (S) versus measured ^4He for both datasets. In 2016, a He sensitivity of 1.1×10^{-12} ccSTP He cps⁻¹ (4.9×10^{-17} mol He cps⁻¹) was obtained. In 2018, the addition of a cryogenic trap and modifications of source parameters allowed us to obtain a better sensitivity of 2.9×10^{-14} ccSTP He cps⁻¹ (1.3×10^{-18} mol He cps⁻¹). The variable n is the number of analyzed fragments.




Figure 6. Evolution of the $^{235}\text{U}/^{238}\text{U}$, $^{230}\text{Th}/^{232}\text{Th}$, and $^{149}\text{Sm}/^{147}\text{Sm}$ ratios obtained for BSP, BSP-ch, and Durango apatite Sp solutions in six different chemistry sessions.

BSP. In addition, the use of Nb capsules (BSP-ch-Nb) also impacts the U and Th budgets and leads to a massive reduction of the ratios (maxima $^{235}\text{U}/^{238}\text{U} = 16.2 \pm 1.0$; $^{230}\text{Th}/^{232}\text{Th} = 4.1 \pm 0.9$ and recorded minima down to $^{235}\text{U}/^{238}\text{U} = 12.0$; $^{230}\text{Th}/^{232}\text{Th} = 1.7$; Fig. 7).

The impact of the niobium capsule on U and Th signals has already been noted by Reiners and Nicolescu (2007) and

reported to be more significant for the Th content. The differences in $^{235}\text{U}/^{238}\text{U}$ and $^{230}\text{Th}/^{232}\text{Th}$ ratios between the different blanks (BSP, BSP-ch, BSP-Ch-Nb) are associated with some isotope fractionation, with a calculated decrease of 4% to 6% for the $^{235}\text{U}/^{238}\text{U}$ ratio and 10 to 20% for the $^{230}\text{Th}/^{232}\text{Th}$ ratio compared to the expected ratios given by the spiked solutions without Nb. The shift for these ratios is systematic but variable from one solution to another, particularly for the $^{230}\text{Th}/^{232}\text{Th}$ ratio, showing values down to 60% of the expected spiked solution ratio. A correction of the impact of the niobium capsule on the U and Th solutions is considered as this effect can result in a shift of the zircon (U–Th) / He ages of up to $\sim 20\%$.

3.4 Durango apatite and FCT zircon (U–Th) / He age reproducibility

The Durango apatite is constantly analyzed in the laboratory to check for the He mass spectrometer sensitivity evolution though time as well as the evolution of (U–Th) / He age and U, Th, and Sm contents. As the dissolution protocol for apatite has a very low impact on the U, Th, and Sm content determination by ICP-MS analysis, the regular measurement of Durango apatite acts as a sensor and allows us to detect any analytical problem. As an example, Fig. 8 presents the values of the (U–Th) / He ages and Th / U and Sm / Th ratios acquired from March to December 2019 by both the quad and VG lines, and values are reported in Table S1 in the Supplement. The mean of the (U–Th) / He age is 31.1 ± 1.4 Ma, in agreement with the age of 31.02 ± 1.01 Ma from McDowell et al. (2005).

A typical mean error of $< 5\%$ (1σ) is obtained for each AHe age by using either the quadrupole or magnetic sec-




Figure 7. Evolution of the $^{235}\text{U}/^{238}\text{U}$ and $^{230}\text{Th}/^{232}\text{Th}$ ratios obtained for BSP, BSP-ch, BSP-ch-Nb, and FCT zircon solutions in different chemistry sessions. 1-Zr, 2-Zr, and 3-Zr refer to the spike solution names used during the dissolution.




Figure 8. Durango (U–Th) / He age dispersion (1σ shown) as a function of Th / U (a) or Sm / Th ratios (b). Durango AHe ages and elemental ratios acquired from March to December 2019 using the quad line (open red circles) or the VG line (filled red circles). Histogram representation of the He ages; Th / U and Sm / Th ratios were constructed using RadialPlotter (Vermeesch, 2009). Th / U and Sm / Th ratios have no unit but are expressed by mass (ng ng^{-1}).




Figure 9. Fish Canyon Tuff zircon (U–Th) / He age dispersion (1σ shown) as a function of the Th / U ratio for the data obtained in 2018 using the quad and VG lines. Histogram representations of the (U–Th) / He ages and Th / U ratios were constructed using RadialPlotter (Vermeesch, 2009). A U / Pb age of 28.5 ± 0.06 Ma has been published by Schmitz and Bowring (2001). Th / U is expressed by mass.

tor mass spectrometers, without any evident difference over a large period of time. This error can be interpreted as the quadratic sum of the errors in the coupled analyses of the U–Th–Sm and He contents and is associated with the calibration. Our (U–Th) / He ages for unknown apatites compare well within error with other laboratories (e.g., Ketcham et al., 2018). In addition, our strategy developed to determine the Ca concentration allows us to obtain the weight of the Durango fragment(s) and thus to calculate the U, Th, and Sm concentrations in parts per million (ppm; Table S1). Mean values of $U = 19 \pm 4$ ppm, $Th = 412 \pm 68$ ppm, and $Sm = 38 \pm 7$ ppm have been obtained, and the U content is similar to that obtained by Schneider et al. (2015) and Yanga et al. (2014).

Fish Canyon Tuff zircon crystals have been analyzed at the GEOPS laboratory as standards for the ZHe method. The U and Th losses during dissolution due to niobium impact are corrected on the determination of the U and Th concentrations in zircons. (U–Th) / He ages were obtained on 57 crystals of FCT zircons analyzed using the quad and VG lines and are reported in Fig. 9 as a function of the Th / U ratio as well as in Table S2.

We obtain a mean age of 27.0 ± 2.6 Ma (1σ) and a mean Th / U ratio of 0.4 on a zircon crystal from the C. W. Naeser collection (K / Ar age of 27.9 ± 0.7 Ma; Naeser et al., 1981) (Fig. 9). The standard dispersion of the ZHe ages is $\sim 9\%$ and is comparable to the natural dispersion observed in the ZHe values given in the literature (e.g., Ault et al., 2019; Guenther et al., 2014; Reiners, 2005). The Th / U dispersion of 37 % also corresponds to the natural dispersion observed in the Th / U ratio of the Fish Canyon zircon standard (e.g., Reiners et al., 2002). The (U–Th) / He age results are comparable with (U–Th) / He literature data that range from 27.3 ± 1.0 to 29.8 ± 2.7 Ma (Dobson et al., 2008; Gleadow et

al., 2015; Reiners et al., 2002; Tagami et al., 2003; Tibari et al., 2016). Th / U ratios vary between 0.42 ± 0.15 (Tagami et al., 2003) and 0.63 ± 0.14 Ma (Tibari et al., 2016) in the literature. The mean (U–Th) / He age obtained in this study is slightly younger, by 5 %, than the U / Pb age of 28.5 ± 0.06 Ma obtained by Schmitz and Bowring (2001), but it is still in good agreement within error bars (Fig. 9). The slight ZHe age difference could be explained by the variability of the measured Fish Canyon Tuff zircon ages as a function of sampling site (Gleadow et al., 2015). A second option is that, since similar ages are obtained by degassing either on the quad or VG line, the slight shift in the (U–Th) / He age may be associated with the He content determination for only 1 % or 2 %, but it is mainly associated with the U and Th content determination and finally with the impact of the niobium precipitation during zircon dissolution. The loss of U, and more specifically of Th, associated with the use of $HF + HNO_3$ (Révillon and Hureau-Mazaudier, 2009; Yokoyama et al., 1999) and of Nb capsules (Reiners and Nicolescu, 2007) has already been taken into account in the blank correction. However, the slightly lower ZHe ages obtained in this study could also be associated with the impact of the zirconium brought into solution that could cause an additional slight loss in Th, which is not considered yet in the blank correction. To estimate the magnitude of this under-correction, additional work should be carried out to fully understand the U and Th isotopic and elemental fractionation during this chemical protocol.

4 Conclusion

This contribution presents the (U–Th) / He analysis protocols developed over the last 10 years at the GEOPS laboratory, Paris-Saclay University, and shares all the empirical

and analytical aspects observed during the different steps of the protocol: sample preparation, mineral hand picking, He analysis, mineral dissolution, and U, Th, and Sm content determination. In the light of our experience, we propose the following:

- a simple method to determine the temperature of the heated metallic (Pt and Nb) capsules that contain the apatite or zircon crystals during laser firing in the range 900–1200 °C using visible light emission wavelengths;
- a method to calibrate He sensitivity using quadrupole and magnetic sector mass spectrometers;
- the protocols to dissolve apatite and zircon crystals and to clean laboratory vessels after chemical digestion;
- the protocol to calibrate the U, Th, and Sm spikes; and
- the method used to track the U, Th, and Sm blank evolution and determine U, Th, and Sm contents.

We adopted the Durango apatite as a standard to perform He calibration and check for He and U–Th–Sm analytical problems; we can thus determine (U–Th) / He ages with an error of less than 5 % (1σ). Our choice is also related to the fact that Durango is an easy-to-use mineral due to its high purity, its rapid dissolution protocol, and the strong reproducibility of its analyses. For the long-term quality control of the (U–Th) / He data, attention needs to be paid to precisely evaluating the drift of blanks through time and that of the (U–Th) / He ages and Th / U ratios (with Sm / Th when possible) obtained on standards (Durango apatite and Fish Canyon Tuff zircon), especially when using quadrupole mass spectrometry.

Appendix A

Table A1. U, Th, Sm, and Ca spike calibration table (ppm and ppb: $\mu\text{g g}^{-1}$ and ng g^{-1}).

Comment	Isotope	Volume	Solutions
Spike solutions preparation			
MR: concentrated mono-elemental mother solution and solid	$^{235}\text{U}_{\text{MR}}$	10 mL	HNO_3 5N 1 ppm (solution)
	$^{230}\text{Th}_{\text{MR}}$	5 mL	HNO_3 5N 0.569 ppm (solution)
	$^{149}\text{Sm}_{\text{MR}}$	100 mL	HNO_3 5N 7 ppm (solution)
	$^{42}\text{Ca}_{\text{MR}}$		2.44 mg $^{42}\text{CaCO}_3$ (solid)
MR1: first-dilution mother spike solution obtained by dilution from the concentrated MR mother solution or solid	$^{230}\text{Th}_{\text{MR1}} = 100$ ppb	10 mL	No need to dilute ^{235}U and ^{149}Sm (already at adapted concentration)
	$^{42}\text{Ca}_{\text{MR1}} = 100$ ppm	10 mL	1.8 mL $^{230}\text{Th}_{\text{MR}}$ + 8.2 mL HNO_3 5N 2.44 mg $^{42}\text{CaCO}_3$ (solid) in 10 mL HNO_3 5N
MR2: second-dilution mother spike solution obtained by mixing–dilution of the concentrated (MR) and first-dilution mother spike (MR1)	$^{235}\text{U}_{\text{MR2}} = 4$ ppb $^{230}\text{Th}_{\text{MR2}} = 4$ ppb $^{149}\text{Sm}_{\text{MR2}} = 4$ ppb $^{42}\text{Ca}_{\text{MR2}} = 1$ ppm	60 mL	240 μL $^{235}\text{U}_{\text{MR}}$ (1 ppm) + 2.4 mL $^{230}\text{Th}_{\text{MR1}}$ (100 ppb) + 34 μL $^{149}\text{Sm}_{\text{MR}}$ (7 ppm) + 598 μL $^{42}\text{Ca}_{\text{MR1}}$ (100 ppm) + 56.727 mL HNO_3 5N
Spiked solutions for spikes calibration			
S: concentrated standard solutions used for spike calibration	$^{238}\text{U}_{\text{S}}$ $^{232}\text{Th}_{\text{S}}$ Sm_{S} natural Ca_{S} natural	125 mL	^{238}U tailored solution HNO_3 5N 1015 ppm ^{232}Th tailored solution HNO_3 5N 993 ppm Natural Sm solution HNO_3 5N 1006 ppm ($^{147}\text{Sm} = 14.99\%$; $^{149}\text{Sm} = 13.82\%$) Natural Ca solution HNO_3 5N 1003 ppm ($^{40}\text{Ca} = 96.94\%$, $^{42}\text{Ca} = 0.65\%$, $^{43}\text{Ca} = 0.14\%$)
U_{S} (III): freshly made multi-step dilutions to obtain 4 ppb ^{238}U mono-elemental standard solution	$^{238}\text{U}_{\text{S}}$ (I)	8 mL	(I) 10 ppm $^{238}\text{U} = 80 \mu\text{L}$ ^{238}U 1015 ppm + 7.92 mL HNO_3 5N
	$^{238}\text{U}_{\text{S}}$ (II)	8 mL	(II) 100 ppb $^{238}\text{U} = 80 \mu\text{L}$ ^{238}U 10 ppm + 7.92 mL HNO_3 5N
	$^{238}\text{U}_{\text{S}}$ (III)	8 mL	(III) 4 ppb $^{238}\text{U} = 320 \mu\text{L}$ ^{238}U 100 ppb + 7.68 mL HNO_3 5N
Th_{S} (III): freshly made multi-step dilutions to obtain 4 ppb ^{232}Th mono-elemental standard solution	$^{232}\text{Th}_{\text{S}}$ (I)	8 mL	(I) 10 ppm $^{232}\text{Th} = 80 \mu\text{L}$ ^{232}Th 993 ppm + 7.92 mL HNO_3 5N
	$^{232}\text{Th}_{\text{S}}$ (II)	8 mL	(II) 100 ppb $^{232}\text{Th} = 80 \mu\text{L}$ ^{232}Th 10 ppm + 7.92 mL HNO_3 5N
	$^{232}\text{Th}_{\text{S}}$ (III)	8 mL	(III) 4 ppb $^{232}\text{Th} = 320 \mu\text{L}$ ^{232}Th 100 ppb + 7.680 mL HNO_3 5N
Sm_{S} (III): freshly made multi-step dilutions to obtain 4 ppb ^{149}Sm mono-elemental standard solution	Sm_{S} (I)	10 mL	(I) 1 ppm Sm nat = 10 μL Sm nat 1006 ppm + 9.99 mL HNO_3 5N
	Sm_{S} (II)	10 mL	(II) 70 ppb Sm nat = 0.7 mL Sm nat 1 ppm + 9.3 mL HNO_3 5N
	Sm_{S} (III)	8 mL	(III) 4 ppb Sm nat = 0.46 mL Sm nat 70 ppb + 7.54 mL HNO_3 5N
SS: freshly spiked standard solutions prepared using MR2 and S solutions	U_{SS} Th_{SS} Sm_{SS} Ca_{SS}	2 mL	50 μL $^{238}\text{U}_{\text{S}}$ (III) (4 ppb) + 50 μL $^{235}\text{U}_{\text{MR2}}$ (4 ppb) + 1.9 mL HNO_3 1N 50 μL $^{232}\text{Th}_{\text{S}}$ (III) (4 ppb) + 50 μL $^{230}\text{Th}_{\text{MR2}}$ (4 ppb) + 1.9 mL HNO_3 1N 50 μL Sm_{S} (III) (4 ppb) + 50 μL $^{149}\text{Sm}_{\text{MR2}}$ (4 ppb) + 1.9 mL HNO_3 1N 5 μL Ca_{S} (1003 ppm) + 50 μL $^{42}\text{Ca}_{\text{MR2}}$ (1 ppm) + 1.945 mL HNO_3 1N

Appendix B: Equations used to determine the ^{238}U , ^{232}Th , ^{147}Sm , and ^{43}Ca contents

The U, Th, Sm, and Ca abundances are deduced from the measured $^{235}\text{U} / ^{238}\text{U}$, $^{230}\text{Th} / ^{232}\text{Th}$, $^{149}\text{Sm} / ^{147}\text{Sm}$, and $^{42}\text{Ca} / ^{43}\text{Ca}$ isotopic ratios and the equations listed below (Eqs. B1–B5) derived from Evans et al. (2005).

$$^{238}\text{U} = \left(m_{\text{spk}}(g) \times \left[^{235}\text{U}_{\text{spk}} \right] \times 10^{-9} \right) \times \frac{1}{\left(\left(\frac{^{235}\text{U}}{^{238}\text{U}} \right)_{\text{spk}} \right)} \times \left(\frac{\left(\frac{^{235}\text{U}}{^{238}\text{U}} \right)_{\text{sp}} - \left(\frac{^{235}\text{U}}{^{238}\text{U}} \right)_{\text{BSP}}}{\left(\frac{^{235}\text{U}}{^{238}\text{U}} \right)_{\text{Nat}} - \left(\frac{^{235}\text{U}}{^{238}\text{U}} \right)_{\text{sp}}} \right) \quad (\text{B1})$$

$$^{232}\text{Th} = \left(m_{\text{spk}}(g) \times \left[^{230}\text{Th}_{\text{spk}} \right] \times 10^{-9} \right) \times \frac{1}{\left(\left(\frac{^{230}\text{Th}}{^{232}\text{Th}} \right)_{\text{spk}} \right)} \times \left(\frac{\left(\frac{^{230}\text{Th}}{^{232}\text{Th}} \right)_{\text{BSP}} - \left(\frac{^{230}\text{Th}}{^{232}\text{Th}} \right)_{\text{sp}}}{\left(\frac{^{230}\text{Th}}{^{232}\text{Th}} \right)_{\text{sp}}} \right) \quad (\text{B2})$$

$$^{147}\text{Sm} = \left(m_{\text{spk}}(g) \times \left[^{149}\text{Sm}_{\text{spk}} \right] \times 10^{-9} \right) \times \frac{1}{\left(\left(\frac{^{149}\text{Sm}}{^{147}\text{Sm}} \right)_{\text{spk}} \right)} \times \left(\frac{\left(\frac{^{149}\text{Sm}}{^{147}\text{Sm}} \right)_{\text{sp}} - \left(\frac{^{149}\text{Sm}}{^{147}\text{Sm}} \right)_{\text{BSP}}}{\left(\frac{^{149}\text{Sm}}{^{147}\text{Sm}} \right)_{\text{Nat}} - \left(\frac{^{149}\text{Sm}}{^{147}\text{Sm}} \right)_{\text{sp}}} \right) \quad (\text{B3})$$

$$^{43}\text{Ca} = \left(m_{\text{spk}}(g) \times \left[^{42}\text{Ca}_{\text{spk}} \right] \times 10^{-9} \right) \times \frac{1}{\left(\left(\frac{^{42}\text{Ca}}{^{43}\text{Ca}} \right)_{\text{spk}} \right)} \times \left(\frac{\left(\frac{^{42}\text{Ca}}{^{43}\text{Ca}} \right)_{\text{sp}} - \left(\frac{^{42}\text{Ca}}{^{43}\text{Ca}} \right)_{\text{BSP}}}{\left(\frac{^{42}\text{Ca}}{^{43}\text{Ca}} \right)_{\text{Nat}} - \left(\frac{^{42}\text{Ca}}{^{43}\text{Ca}} \right)_{\text{sp}}} \right) \quad (\text{B4})$$

The following abbreviations are applied: spk – spike, sp – sample, BSP – spiked blank, Nat – natural ratio, i.e., $^{235}\text{U} / ^{238}\text{U} = 0.00725268$, $^{149}\text{Sm} / ^{147}\text{Sm} = 0.92$ ($^{147}\text{Sm} = 14.99\%$ and $^{149}\text{Sm} = 13.8\%$), and $^{42}\text{Ca} / ^{43}\text{Ca} = 4.79$ ($^{42}\text{Ca} = 0.647\%$ and $^{43}\text{Ca} = 0.135\%$). We obtain abundances for each sample in nanograms (10^{-9} g), with spike concentrations in parts per billion (ppb). In addition, the same equations can be reversed to determine the concentration of the spike isotopes ^{235}U , ^{230}Th , ^{149}Sm , and ^{42}Ca .

Guenther et al. (2016) have already performed complete work on determining apatite and zircon weights by measuring the Ca and Zr contents. Here, we obtain the weight of the apatite grain(s) from the measurement of ^{43}Ca (ng) and the composition of a pure fluoroapatite ($\text{Ca}_5(\text{PO}_4)_3\text{F}$) containing 40 wt % Ca in one apatite crystal. Thus, the apatite weight is given by Eq. (B5).

$$\text{weight } (\mu\text{g}) = \frac{^{43}\text{Ca} \times 10^{-9}}{0.135/100} \times \frac{1}{0.4 \times 10^{-6}} \quad (\text{B5})$$

The factor 0.135 refers to the natural isotope abundance of ^{43}Ca .

Appendix C: Equations used for the (U–Th) / He age reduction

The (U–Th) / He age (Ma) is calculated assuming a linear production of ^4He with time using the determined U, Th, and Sm abundances and Eqs. (C1) and (C2).

$$\frac{(\text{U–Th})}{\text{He}} \text{age (Ma)} = \frac{^4\text{He}}{P^* \times 10^6} \quad (\text{C1})$$

P^* is the instantaneous production of ^4He , i.e., the amount of ^4He produced in 1 year, and ^4He is the measured ^4He abundance (both in ccSTP). P^* is calculated using the following equation.

$$P^* = \left(\frac{8 \times [^{238}\text{U}] \times 10^{-9}}{238} \times \lambda_{^{238}\text{U}} + \frac{7 \times [^{238}\text{U}] \times 10^{-9}}{235 \times 137.88} \times \lambda_{^{235}\text{U}} + \frac{6 \times [^{232}\text{Th}] \times 10^{-9}}{232} \times \lambda_{^{232}\text{Th}} + \frac{[^{147}\text{Sm}] \times 10^{-9}}{147} \times \lambda_{^{147}\text{Sm}} \right) \times 22414 \quad (\text{C2})$$

Here, $\lambda_{^{238}\text{U}}(\text{yr}^{-1}) = \text{Ln}(2)/(4.47 \times 10^9)$, $\lambda_{^{235}\text{U}}(\text{yr}^{-1}) = \text{Ln}(2)/(7.04 \times 10^8)$, $\lambda_{^{232}\text{Th}}(\text{yr}^{-1}) = \text{Ln}(2)/(1.40 \times 10^{10})$, and $\lambda_{^{147}\text{Sm}}(\text{yr}^{-1}) = \text{Ln}(2)/(1.06 \times 10^{11})$, where $[^{238}\text{U}]$, $[^{232}\text{Th}]$, and $[^{147}\text{Sm}]$ are the measured abundances in nanograms. The value of 137.88 is the $^{238}\text{U} / ^{235}\text{U}$ natural isotopic ratio of uranium, and 1 mole occupies 22414 cm^3 at the standard pressure and temperature conditions of 273.15 K and 1 atm (ccSTP). While ccSTP is not the SI unit for amount of substance, it is historically used by the noble gas community, and conversion to moles is easy using the conversion factor of $1 \text{ mol} = 22414 \text{ ccSTP}$.

Code and data availability. The Qt_LFT software and an Excel[®] file example are available in the Supplement, as are Tables S1 and S2 that present Durango and Fish Canyon Tuff zircon (U–Th) / He data. The data reduction Excel[®] workbook is available upon request to Alexis Derycke (alexis.derycke@hotmail.com).

Supplement. The supplement related to this article is available online at: <https://doi.org/10.5194/gchron-3-351-2021-supplement>.

Author contributions. CG, RPJ, PS, LTG, JM, and DB designed the experiments. CG, RPJ, FA, AD, CS, FH, GM, and GD participated in the data acquisition. LTG developed the Monte Carlo simulation Qt_LFT software and AD the Excel[®] workbook automatization software. CG prepared the paper with contributions from all co-authors.

Competing interests. The authors declare that they have a possible conflict of interest as Cécile Gautheron is a member of the editorial board of the journal.

Acknowledgements. The analytical work and the He and VG line building have been funded thanks to the INSU-Relief, FORPRO, and Tellus programs, as well as INSU mi-lourd, Division de la recherche of the Paris Sud University, and the ERM Paris Sud program under grants ANR-06-JCJC-0079 and ANR-12-NS06-0005-01 (HeDiff projects). Maurice Pagel is warmly thanked for his help in funding the cryogenic trap. IUT-Mesures-physiques internship trainees Pierre Marty, Stephane Lemaire, Etienne Cornier, Moussa Di Giacomo, Guillaume Ya, Pierre Boutteville, Maxime Form, Nathasha Etienne, Claire Morelière, Hugo Vicente, and Benoit Canguilhem are warmly thanked for their work in LabView programming of the different parts of the quad and VG lines. We thank Peter Reiners and Uttam Chowdhury for sharing zircon fractions and their knowledge of zircon dissolution protocols, and we thank James R. Metcalf for the location of the FCT outcrop and advice. We warmly thank Peter C. Hackspacher for the donation of gem-quality Durango apatite crystals and Charles W. Naeser for the donation of Fish Canyon Tuff zircon fractions. Manuel Moreira is thanked for his help in building the valve-controlling electronics and building the VG line. Jean Louis Birck is thanked for the loan of the external filament extinction pyrometer. Pierre Burckel, Eric Douville, and Louise Bordier are thanked for the U–Th–Sm analyses at IGP and LSCE. The anonymous reviewers are thanked for their constructive reviews, and the associate editor Marissa Tremblay and editor Greg Balco are thanked for review handling.

Financial support. This research has been supported by the Agence de la recherche (HeDiff project, grant nos. ANR-06-JCJC-0079 and ANR-12-NS06-0005-01), INSU-Relief, FORPRO, and Tellus programs, as well as INSU mi-lourd, Division de la recherche of the Paris Sud University, and the ERM Paris Sud program.

Review statement. This paper was edited by Marissa Tremblay and reviewed by two anonymous referees.

References

- Allard, T., Gautheron, C., Bressan-Riffel, S., Balan, E., Selo, M., Fernandes, B. S., Pinna-Jamme, R., Derycke, A., Morin, G., Taitson Bueno, G., and Do Nascimento, N. R.: Combined dating of goethites and kaolinites from ferruginous duricrusts. Deciphering the Late Neogene erosion history of Central Amazonia, *Chem. Geol.*, 479, 136–150, <https://doi.org/10.1016/j.chemgeo.2018.01.004>, 2018.
- Ault, A. K., Gautheron, C., and King, G. E.: Innovations in (U–Th) / He, fission-track, and trapped-charge thermochronometry with applications to earthquakes, weathering, surface-mantle connections, and the growth and decay of mountains, *Tectonics*, 38, 3705–3739, <https://doi.org/10.1029/2018TC005312>, 2019.
- Brown, R. W., Beucher, R., Roper, S., Persano, C., Stuart, F., and Fitzgerald, P.: Natural age dispersion arising from the analysis of broken crystals, Part I. Theoretical basis and implications for the apatite (U–Th) / He thermochronometer, *Geochim. Cosmochim. Acta*, 122, 478–497, <https://doi.org/10.1016/j.gca.2013.05.041>, 2013.
- Burnard, P. G. and Farley, K. A.: Calibration of pressure-dependent sensitivity and discrimination in Nier-type noble gas ion sources, *Geochem. Geophys. Geosy.*, 1, 1022, <https://doi.org/10.1029/2000GC000038>, 2000.
- Cooperdock, E. H. G., Ketcham, R. A., and Stockli, D. F.: Resolving the effects of 2-D versus 3-D grain measurements on apatite (U–Th) / He age data and reproducibility, *Geochronology*, 1, 17–41, <https://doi.org/10.5194/gchron-1-17-2019>, 2019.
- Cros, A., Gautheron, C., Pagel, M., Berthet, P., Tassan-Got, L., Douville, E., Pinna-Jamme, R., and Sarda, P.: ⁴He behavior in calcite filling viewed by (U–Th) / He dating, ⁴He diffusion and crystallographic studies, *Geochim. Cosmochim. Acta*, 125, 414–432, <https://doi.org/10.1016/j.gca.2013.09.038>, 2014.
- Dobson, K. J., Stuart, F. M., and Dempster, T. J.: U and Th zonation in Fish Canyon Tuff zircons: Implications for a zircon (U–Th) / He standard, *Geochim. Cosmochim. Acta*, 72, 4745–4755, <https://doi.org/10.1016/j.gca.2008.07.015>, 2008.
- Evans, N. J., Byrne, J. P., Keegan, J. T., and Dotter, L. E.: Determination of uranium and thorium in zircon, apatite, and fluorite: Application to laser (U–Th) / He thermochronology, *J. Ana. Chem.*, 60, 1300–1307, 2005.
- Farley, K. A.: Helium diffusion from apatite: general behavior as illustrated by Durango fluorapatite, *J. Geophys. Res.*, 105, 2903–2914, 2000.
- Farley, K. A.: (U–Th) / He dating: Techniques, calibrations, and applications, in: *geochemistry*, edited by: Porcelli, D., Ballentine, C. J., and Wieler, R., Noble Gases in Geochemistry and Cosmochemistry, 819–844, <https://doi.org/10.2138/rmg.2002.47.18>, 2002.
- Farley, K. A., Wolf, R. A., and Silver, L. T.: The effects of long alpha-stopping on (U–Th)/He ages, *Geochim. Cosmochim. Acta*, 21, 4223–4229, [https://doi.org/10.1016/S0016-7037\(96\)00193-7](https://doi.org/10.1016/S0016-7037(96)00193-7), 1996.
- Fechtig, H. and Kalbitzer, S.: Diffusion of argon in potassium-bearing solids, in: *Potassium-Argon Dating*, edited by:

- Zahringer, O. S. a. J., Springer, Berlin, Heidelberg, Germany, 68–106, 1966.
- Flowers, R., Ketcham, R. A., Shuster, D., and Farley, K. A.: Apatite (U–Th) / He thermochronology using a radiation damage accumulation and annealing model, *Geochim. Cosmochim. Acta*, 73, 2347–2365, <https://doi.org/10.1016/j.gca.2009.01.015>, 2009.
- Foeken, J. P. T., Stuart, F. M., Dobson, K. J., Persano, C., and Vilbert, D.: A diode laser system for heating minerals for (U–Th) / He chronometry, *Geochem. Geophys.*, 7, 1–9, <https://doi.org/10.1029/2005GC001190>, 2006.
- Gautheron, C. and Tassan-Got, L.: A Monte Carlo approach of diffusion applied to noble gas/ helium thermochronology, *Chem. Geol.*, 273, 212–224, <https://doi.org/10.1016/j.chemgeo.2010.02.023>, 2010.
- Gautheron, C. and Zeitler, P. K.: Noble Gases Deliver Cool Dates from Hot Rocks, *Elements*, 16, 303–309, <https://doi.org/10.2138/gselements.16.5.303>, 2020.
- Gautheron, C., Tassan-Got, L., Ketcham, R. A., and Dobson, K. J.: Accounting for long alpha-particle stopping distances in (U–Th–Sm) / He geochronology: 3D modeling of diffusion, zoning, implantation, and abrasion, *Geochim. Cosmochim. Acta*, 96, 44–56, <https://doi.org/10.1016/j.gca.2012.08.016>, 2012.
- Gautheron, C., Djimbi, D. M., Roques, J., Balout, H., Ketcham, R. A., Simoni, E., Pik, R., Seydoux-Guillaume, A. M., and Tassan-Got, L.: A multi-method, multi-scale theoretical study of He and Ne diffusion in zircon, *Geochim. Cosmoch. Acta*, 268, 348–367, <https://doi.org/10.1016/j.gca.2019.10.007>, 2020.
- Gautheron, C. E., Tassan-Got, L., Barbarand, J., and Pagel, M.: Effect of alpha-damage annealing on apatite (U–Th) / He thermochronology, *Chem. Geol.*, 266, 166–179, <https://doi.org/10.1016/j.chemgeo.2009.06.001>, 2009.
- Gerin, C., Gautheron, C., Oliviero, E., Bachelet, C., Djimbi, D. M., Seydoux-Guillaume, A. M., Tassan-Got, L., Sarda, P., Roques, J., and Garrido, F.: Influence of vacancy damage on He diffusion in apatite investigated at atomic to mineralogical scales, *Geochim. Cosmoch. Acta*, 197, 87–103, <https://doi.org/10.1016/j.gca.2016.10.018>, 2017.
- Gleadow, A., Harrison, M., Kohn, B. P., Lugo-Zazuta, R., and Phillips, D.: The Fish Canyon Tuff: A new look at an old low-temperature thermochronology standard, *Earth Planet. Sci. Lett.*, 424, 95–108, <https://doi.org/10.1016/j.epsl.2015.05.003>, 2015.
- Glotzbach, C., Lang, K. A., Avdievitch, N. N., and Ehlers, T.: Increasing the accuracy of (U–Th–Sm) / He dating with 3D grain modelling, *Chem. Geol.*, 506, 113–125, <https://doi.org/10.1016/j.chemgeo.2018.12.032>, 2019.
- Goldsmith, A. S., Ketcham, R. A., and Stockli, D.: Simulating effects of heterogeneous ⁴He concentration profiles and radiation damage annealing on whole-grain zircon diffusivity analyses, *Geochim. Cosmochim. Acta*, 284, 239–253, <https://doi.org/10.1016/j.gca.2020.05.033>, 2020.
- Guenther, W., Reiners, P. W., Ketcham, R., Nasdala, L., and Gierster, G.: Helium diffusion in natural zircon: radiation damage, anisotropy, and the interpretation of zircon (U–Th) / He thermochronology, *Am. J. Sci.*, 313, 145–198, <https://doi.org/10.2475/03.2013.01>, 2013.
- Guenther, W. R., Reiners, P. W., DeCelles, P. G., and Kendall, J.: Sevier belt exhumation in central Utah constrained from complex zircon (U–Th) / He data sets: Radiation damage and He inheritance effects on partially reset detrital zircons, *GSA Bulletin*, 1, B31032, <https://doi.org/10.1130/B31032.1>, 2014.
- Guenther, W. R., Reiners, P. W., and Chowdhury, U.: Isotope dilution analysis of Ca and Zr in apatite and zircon (U–Th) / He chronometry, *Geochem. Geophys.*, 17, 1623–1640, <https://doi.org/10.2475/03.2013.01>, 2016.
- Herman, F., Braun, J., Senden, T. J., and Dunlap, W. J.: (U–Th) / He chronometry: Mapping 3D geometry using micro-X-ray tomography and solving the associated production-diffusion equation, *Chem. Geol.*, 242, 126–136, <https://doi.org/10.1016/j.chemgeo.2007.03.009>, 2007.
- House, M. A., Farley, K. A., and Stockli, D.: Helium chronometry of apatite and titanite Nd-YAG laser heating, *Earth Planet. Sci. Lett.*, 183, 365–368, 2000.
- Ketcham, R. A., Gautheron, C., and Tassan-Got, L.: Accounting for long alpha-particle stopping distances in (U–Th–Sm) / He geochronology: refinement of the baseline case, *Geochim. Cosmochim. Acta*, 75, 7779–7791, <https://doi.org/10.1016/j.gca.2011.10.011>, 2011.
- Ketcham, R. A., van der Beek, P., Barbarand, J., Bernet, M., and Gautheron, C.: Reproducibility of Thermal History Reconstruction From Apatite Fission-Track and (U–Th) / He Data, *Geochem. Geophys.*, 19, 2411–2436, <https://doi.org/10.1029/2018GC007555>, 2018.
- McDowell, F. W., McIntosh, W. C., and Farley, K. A.: A precise ⁴⁰Ar–³⁹Ar reference age for the Durango apatite (U–Th) / He and fission-track dating standard, *Chem. Geol.*, 214, 249–263, <https://doi.org/10.1016/j.chemgeo.2004.10.002>, 2005.
- Naeser, C. W., Zimmermann, R. A., and Cebula, G. T.: Fission-track dating of apatite and zircon: An interlaboratory comparison, *Nucl. Tracks*, 5, 65–72, [https://doi.org/10.1016/0191-278X\(81\)90027-5](https://doi.org/10.1016/0191-278X(81)90027-5), 1981.
- Reiners, P. W.: Zircon (U–Th) / He thermochronometry, *Reviews in Mineralogy and Geochemistry*, 58, 151–179, <https://doi.org/10.2138/rmg.2005.58.6>, 2005.
- Reiners, P. W. and Brandon, M. T.: Using thermochronology to understand orogenic erosion, *Ann. Rev. Earth. Planet. Sci.*, 34, 419–466, <https://doi.org/10.1146/annurev.earth.34.031405.125202>, 2006.
- Reiners, P. W. and Farley, K. A.: Influence of crystal size on apatite (U + Th) / He thermochronology: an example from the Bighorn Mountains, Wyoming, *Earth Planet. Sci. Lett.*, 188, 413–420, [https://doi.org/10.1016/S0012-821X\(01\)00341-7](https://doi.org/10.1016/S0012-821X(01)00341-7), 2001.
- Reiners, P. W. and Nicolescu, S.: Measurement of parent nuclides for (U–Th) / He chronometry by solution sector ICP-MS, ARHDL Report 1, 1–33, available at: https://www.geo.arizona.edu/~reiners/arhdl/arhdl_reports.htm (last access: 31 May 2021), 2007.
- Reiners, P. W., Farley, K. A., and Hiskes, H. J.: He diffusion and (U–Th) / He thermochronology of zircon: initial results from Fish Canyon Tuff and Gold Butte, *Tectonophysics*, 349, 297–308, [https://doi.org/10.1016/S0040-1951\(02\)00058-6](https://doi.org/10.1016/S0040-1951(02)00058-6), 2002.
- Réveillon, S. and Hureau-Mazaudier, D.: Improvements in Digestion Protocols for Trace Element and Isotope Determinations in Stream and Lake Sediment Reference Materials (JSd-1, JSd-2, JSd-3, Jlk-1 and LKSD-1), *Geostand. Geoanal. Res.*, 33, 397–413, <https://doi.org/10.1111/j.1751-908X.2009.00008.x>, 2009.
- Schmitz, M. D. and Bowring, S. A.: U–Pb zircon and titanite systematics of the Fish Canyon Tuff: an assessment of high

- precision U–Pb geochronology and its application to young volcanic rocks, *Geochim. Cosmochim. Acta*, 65, 2571–2587, [https://doi.org/10.1016/S0016-7037\(01\)00616-0](https://doi.org/10.1016/S0016-7037(01)00616-0), 2001.
- Schneider, S., Hammerschmidt, K., Rosenberg, C. L., Gardes, A., Frei, D., and Bertrand, A.: U–Pb ages of apatite in the western Tauern Window (Eastern Alps): Tracing the onset of collision-related exhumation in the European plate, *Earth Planet. Sci. Lett.*, 418, 53–65, <https://doi.org/10.1016/j.epsl.2015.02.020>, 2015.
- Shuster, D., Flowers, R., and Farley, K. A.: The influence of natural radiation damage on helium diffusion kinetics in apatite, *Earth Planet. Sci. Lett.*, 249, 148–161, <https://doi.org/10.1016/j.epsl.2006.07.028>, 2006.
- Tagami, T., Farley, K. A., and Stockli, D. F.: (U + Th) / He geochronology of single zircon grains of known Tertiary eruption age, *Earth Planet. Sci. Lett.*, 207, 57–67, [https://doi.org/10.1016/S0012-821X\(02\)01144-5](https://doi.org/10.1016/S0012-821X(02)01144-5), 2003.
- Tibari, B., Tibari, B., Vacherat, A., Stab, M., Pik, R., Yeghicheyan, D., and Hild, P.: An Alternative Protocol for Single Zircon Dissolution with Application to (U–Th–Sm) / He Thermochronometry, *Geostand. Geoanal. Res.*, 40, 365–375, <https://doi.org/10.1111/j.1751-908X.2016.00375.x>, 2016.
- Vermeesch, P.: RadialPlotter: A Java application for fission track, luminescence and other radial plots, *Radiat. Meas.*, 44, 409–410, 2009.
- Yanga, Y.-H., Wua, F.-Y., Yanga, J.-H., Chew, D.M., Xie, L.-W., Chua, Z.-Y., Zhang, Y.-B., and Huang, C.: Sr and Nd isotopic compositions of apatite reference materials used in U–Th–Pb geochronology, *Chem. Geol.*, 385, 35–55, <https://doi.org/10.1016/j.chemgeo.2014.07.012>, 2014.
- Yokoyama, T., Makishima, A., and Nakamura, E.: Evaluation of the coprecipitation of incompatible trace elements with fluoride during silicate rock dissolution by acid digestion, *Chem. Geol.*, 157, 175–187, [https://doi.org/10.1016/S0009-2541\(98\)00206-X](https://doi.org/10.1016/S0009-2541(98)00206-X), 1999.
- Ziegler, J. F.: SRIM-2008 The stopping range of ions in matter, United States Naval Academy, Annapolis, USA, 2008.